

# Investigating Your Watershed: Detection of Fecal Coliform Bacteria

*Karin Readel*

Interdisciplinary Science Program  
University of Maryland, Baltimore County  
111 Physics Bldg.  
1000 Hilltop Circle  
Baltimore, MD 21250  
Phone: 410-455-1425  
Fax: 410-455-1072  
readel@umbc.edu

Karin Readel is a lecturer and the coordinator of the Interdisciplinary Science Program at the University of Maryland, Baltimore County. She is developing and teaching a new laboratory-based science course for non-majors entitled “Water: An Interdisciplinary Study.”

**Reprinted From:** Readel, K. 2000. Investigating your watershed: Detection of fecal coliform bacteria. Pages 400-401, *in* Tested studies for laboratory teaching, Volume 22 (S. J. Karcher, Editor). Proceedings of the 22nd Workshop/Conference of the Association for Biology Laboratory Education (ABLE), 489 pages.

- Copyright policy: <http://www.zoo.utoronto.ca/able/volumes/copyright.htm>

Although the laboratory exercises in ABLE proceedings volumes have been tested and due consideration has been given to safety, individuals performing these exercises must assume all responsibility for risk. The Association for Biology Laboratory Education (ABLE) disclaims any liability with regards to safety in connection with the use of the exercises in its proceedings volumes.

©2001 Karin Readel

Fecal coliform bacteria, members of the family Enterobacteriaceae, are found in the feces of warm-blooded animals, including humans. Though not harmful themselves, large numbers of fecal coliform bacteria have been associated with the presence of other pathogenic organisms. These disease-causing organisms are often very hard to detect due to their small numbers, and low survival rate outside their host organisms. Fecal coliform bacteria, on the other hand, can be detected very easily, and are therefore useful indicators of contaminated water.

Micrology Laboratories has developed a simple method for fecal coliform detection called “Coliscan EasyGel®”. Each testing unit consists of a specially treated Petri dish and a bottle of liquid medium. A small sample of the water to be tested is added to the media and poured into the Petri dish. Ions in the Petri dish coating react with the liquid causing it to gel

completely within about 30-40 minutes. Plates are incubated from 24-48 hours (temperature is not critical) at which time colonies can be counted. The media contains two dye- linked sugars that are activated through bacterial enzymes. The enzyme  $\beta$ -galactosidase is produced by general coliform bacteria, including *E. coli* (a fecal coliform), and it results in a pink colony color. The enzyme  $\beta$ -glucuronidase is only produced by fecal coliforms and results in a blue-green color. Therefore fecal coliform colonies will appear as a dark purplish blue as a result of production of both of these dyes, and can easily be distinguished from the pink general coliform colonies.

This assay is one of many exercises used in a non-major course at the University of Maryland, Baltimore County called *Water: An Interdisciplinary Study*. Students work in groups throughout the semester, performing various lab exercises dealing with the chemical, physical, and biological properties of water. They then use these skills to design and complete a group project based on water quality. Many students choose to investigate the water quality of one of the ponds or streams on campus. The detection of high levels of fecal coliform bacteria in one of these ponds, especially after heavy rains, prompted a lively discussion in lecture. It also provided an opportunity for students to look at the drainage patterns within the local watershed and to suggest possible sources of the pollution.

**Ordering Information:** Micrology Laboratories  
206 W. Lincoln Ave.  
Goshen, IN 46526-3219  
Phone: 1-888-327-9435  
Fax: 1-219-533-3370  
Email: [micrologylabs@juno.com](mailto:micrologylabs@juno.com)  
Web: <http://www.micrologylabs.com>