# The Red and White Yeast Lab: An Introduction to Science As A Process

## **Brian White**

Biology Department University of Massachusetts Boston 100 Morrissey Blvd Boston MA 02125-3393 617-287-6630 brian.white@umb.edu

Brian White teaches General Biology I and II at University of Massachusetts, Boston; he gives the lectures, trains the TAs, and develops the lab exercises. His research interests include development and evaluation of educational multimedia and lab exercises.

**Reprinted From:** White, B. 2000. The red and white yeast lab: An introduction to science as a process. Pages 451-452, *in* Tested studies for laboratory teaching, Volume 21 (S. J. Karcher, Editor). Proceedings of the 21st Workshop/Conference of the Association for Biology Laboratory Education (ABLE), 509 pages.

- Copyright policy: http://www.zoo.utoronto.ca/able/volumes/copyright.htm

Although the laboratory exercises in ABLE proceedings volumes have been tested and due consideration has been given to safety, individuals performing these exercises must assume all responsibility for risk. The Association for Biology Laboratory Education (ABLE) disclaims any liability with regards to safety in connection with the use of the exercises in its proceedings volumes.

Red and White Yeast lab is an experimental system designed for investigation by introductory level biology and non-biology students. The experimental question and techniques have been readily mastered by students in the 6th grade and older. This lab activity focuses on the process of scientific inquiry rather than a particular answer. The teacher's role is that of a moderator; the students design, carry out, and interpret their experiments themselves. The most crucial part of the experience is the classroom discussions which bring out the important features of the interplay between theory and experiment.

The experimental system is an engineered strain of bakers' yeast. When a large patch of these cells is grown on a Petri plate, the patch develops a red center and a white edge. The students' task is to find out what they can about why the center is red and the edge is white. Their tools are fresh plates and sterile toothpicks. Although the system and tools are simple, they allow a wide range of possible models, experiments, and pathways to results.

The students begin by brainstorming about what possible mechanisms could explain this phenomenon. They then do a first round of experiments. At the second meeting, they get their

### Volume 21: Mini Workshops

plates and discuss the results - refining their questions as they see fit. They then do a final round of experiments. At a third and final meeting, they discuss both rounds of results and attempt to reach a consensus about what they can conclude.

#### Schedule

Week 0: In lab: • Look at pictures & video on web site & make hypotheses

Week 1: Pre-lab: • One hypothesis

In lab: • Introduction to yeast

• Look at patches on plates

• Discuss models in groups, then as class

• Introduction to techniques

• Discuss experiments as groups, then as class

• Do Round I of experiments

Week 2: In lab: • Look at results in groups

• Discuss results in groups, then as class

• Discuss Round II experiments as groups, then as class

• Do Round II of experiments

Week 3: In lab: • Look at results in groups

• Discuss results in groups, then as class

Lab Report: • Explain clearly an interaction between hypothesis and

experiment

#### Additional resources:

- The yeast strain is available from the ATCC (www.atcc.org), Accession # 201926 and will soon be available from Carolina Biological Supply Co. (www.carolina.com).
- This lab is part of the ResearchLink 2000 Project (www.researchlink.ferris.edu).
- Pictures of experimental results can be found at: http://omega.cc.umb.edu/~bwhite/yeast/yeast.html
- Detailed descriptions of the experiments, methods, and other resources can be found on the web (http://intro.bio.umb.edu/RL2000/R\_W\_Yeast/).

#### Reference

White, B.T. (1999). The red and white yeast lab. The American Biology Teacher, 61 (8), 600-604.